Simply Practical.

The Hill-Rom[®] CareAssist[®] ES Hospital Bed



Enhancing outcomes for patients and their caregivers:



CareAssist[®] ES Bed: Simply practical.

The CareAssist® ES bed is our foundational Med-Surg bed. We've designed it especially for busy caregivers who need to accomplish a multitude of tasks in a minimum amount of time. A number of one-and-done features distinguish the CareAssist[®] ES bed as one of the most practical beds available.

Do more with easy, one-step features.

- + Bedside controls offer the Dining Chair[®] one-button feature
- + Four-corner brake and steer for easy mobility
- + Built-in bed extender lets a single caregiver add length to the bed
- + Line-of-Site® angle indicator lets nurses monitor a patient's head angle at a glance

Protect your patients from skin breakdown.

Led by experienced clinical professionals, our comprehensive Safe Skin[®] program for the prevention and treatment of skin breakdown helps improve outcomes that are designed to be repeatable and sustainable.

The CareAssist[®] ES bed contributes toward this effort with its SlideGuard[®] feature. which helps protect your patient's skin during repositioning by keeping his or her hips in place while the back of the bed elevates.

The bed's open architecture allows you to use almost any commercially available surface, offering maximum flexibility to care for patients at every level of risk for skin breakdown.

The centrally located fifth wheel adds maneuverability to let you steer the bed easily during transport.

而

WE DESIGNED THE CAREASSIST® ES BED ESPECIALLY FOR **BUSY CAREGIVERS** WHO NEED TO ACCOMPLISH A MULTITUDE **OF TASKS** IN A MINIMUM AMOUNT OF TIME.



CareAssist[®] ES Bed: Simply Practical.

At Hill-Rom, we understand that Med-Surg caregivers must possess a multitude of capabilities to treat patients effectively. The practical CareAssist® ES bed and its easy, one-step features is especially designed for busy caregivers who need to accomplish a multitude of tasks in a minimum amount of time.





One-button Dining Chair® position enables caregivers to place patient in the chair position while the bed remains at a low height.



Point-of-Care[®] siderail controls are conveniently located in a central location, on both sides of the bed, where patient and caregiver interaction occurs.



Built-in Bed Extender allows caregivers to easily add length to the bed.

CareAssist[®] ES Bed: Simply Practical.

The CareAssist[®] ES bed is designed to make bedside tasks easy and fast. Caregivers can enjoy the benefits of a bed that delivers essential and dependable technologies to meet the fundamental needs of the Med-Surg care setting.



Dual-locking brake casters with an audible "brake not set" alarm that reminds caregivers to secure the bed.




With the battery backup, there's no need to find an outlet when repositioning the bed (high/low position, head-of-bed, knee section and foot angle) during patient transfer.



SlideGuard[®] bed articulation technology helps support patient positioning.







SafeView[®] Alert System keeps staff informed of a patient's bed status.





The CareAssist® ES bed connects to your hospital system to send data to NaviCare[®] Nurse Call, EMRs or other clinical applications.

CareAssist[®] ES Bed: Simply Practical.

			Frame length with head/footboard (w/ bumpers)	90" (228.60 cm)
IntelliDrive® Powered Transport*	Line-of-Site [®] Feature		Sleep deck	
	Built-in Bed Extender		width normal length	36" (91.44 cm) 80" (203.20 cm)
	Open Architecture (Range of available mattr SlideGuard® Technology Frame	esses)	extended length	84" (213.36 cm)
	Ability to send data to Nurse Call, EMR or oth	er Clinical Applications	Bed clearance under frame	7.5" (19.05 cm)
			Optional scale system	
Chair Position			accuracy	+/- 2.2 lbs or 1.1% of patient weight w/e
	A			is greater
			weigh in any position scale repeatability	yes 0.2% of patient weight
	1			0.270 of patient weight
Integrated Scale*				
		(The second sec		
Battery Backup		and the second s		
				and the second se
S				
Three-Mode Bed Exit System				and a second
(Patient Position Monitor)*			and the second division of the second divisio	1
		· · · · · · · · · · · · · · · · · · ·		
E			Prost	
Emergency CPR ·····				
	All Contractions			
				1
Point-of-Care [®] ····································	A PA			1 - C
Siderail Controls				
		- Con		
			-	-
		A	A CONTRACT	~ ~
OneStep [®] Operation Siderails				843
				UL
*Optional Features				

Technical Specifications

Frame width

40" (102 cm)

Low position floor to deck	15.75" (40 cm)	
High position floor to deck	32.5" (82.55 cm)	
Max head elevation	65°	
Max Trend and Reverse Trend	16°	
Caster size	6" (15 cm)	
Weight limit (Safe working load)	500 lbs	



The pursuit of greater quality never ends.

Hill-Rom continuously innovates to ensure that we're providing the right tools that help you deliver the highest quality care. We take enormous pride in constructing equipment that enhances the caregiver experience, and we're always listening to our customers to see how we can build greater simplicity and ease of use into our products. It's just one indication of our commitment to the future and will drive even higher quality standards across our company.

To experience the difference, contact your Hill-Rom representative at **800-445-3730** or visit **www.hill-rom.com**.



Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products. ©2014 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

ORDER NUMBER
CTG593ra rev 9
09-MAY-2014
ENG – US

Zinc Number US-AC-1013-0499(1)

</t

www.hill-rom.com

USA 800-445-3730 CANADA 800-267-2337 Enhancing outcomes for patients and their caregivers:

